


# Umbilical Tape

Sterile Nonabsorbable Surgical Cotton



- Soft white cotton
- Consistency you can depend on
- Economical
- Available in 1 dozen vs. 3 dozen

Product	Catalog No.	Length	Per Packet	Per Box
Umbilical Tape	CP10A	1/8" x 18"	2 Strands	1 Dozen
Umbilical Tape	CP11A	1/8" x 30"	2 Strands	1 Dozen
Umbilical Tape	CP12A	1/8" x 36"	3 Strands	1 Dozen

# Bone Wax

CP Medical's Bone Wax is:

- Softer Formulation
- Easy to form and to apply
- 100% USP/NF sterile components
- All natural softener
- Easy to handle



Bone Wax is used to control bleeding on bone surfaces in:

- Thoracic Surgery
- Neurosurgery
- Orthopedic Surgery
- Oral & Maxillofacial Surgeries

Product	Catalog No.	Packaged
Bone Wax	CPB31A	1 Dozen