



CP-FIBER

Ultra high molecular weight polyethylene braided non-absorbable suture



Available with Blue or Black Tracer for easier identification, sizing, motion detection, and arthroscopic visualization

Polyethylene Braided
Polyester Core

STRENGTH

- 2 times stronger than polyester size for size
- 15 times Stronger than Steel on a weighted basis

VARIETY

Available in large variety of suture sizes:
4/0, 3/0, 2/0, 1, 2, 3, 4, 5.

Also available, with or without needles, sterile or non-sterile, and multiple color tracer options.

DIAMETER

Unlike other comparable products, CP-Fiber meets USP Standards for diameter

CP Medical
Portland, OR USA



STERILE PRODUCT

PRODUCT CODE	SIZE	NEEDLE	COLOR	TRACER	LENGTH	PACKAGED
CPF-02	2	NM0-6	White	•	36"	1 Dozen
CPF-02C	2	NM0-6	White	Blue	36"	1 Dozen
CPF-02K	2	NM0-6	White	Black	36"	1 Dozen
CPF-05	5	•	White	•	36"	1 Dozen
CPF-05-SA	5	NCCS	White	•	36"	1 Dozen
BCPF-02-25	2	•	White	•	36"	Bulk - 25 pack
BCPF-02C-25	2	•	White	Blue	36"	Bulk - 25 pack
BCPF-02K-25	2	•	White	Black	36"	Bulk - 25 pack
BCPF-20-25	2/0	•	White	•	36"	Bulk - 25 pack

NON-STERILE PRODUCT

PRODUCT CODE	SIZE	NEEDLE	COLOR	TRACER	LENGTH	PACKAGED
CPF-02-NS	2	NM0-6	White	•	36"	1 Dozen
CPF-02C-NS	2	NM0-6	White	Blue	36"	1 Dozen
CPF-02K-NS	2	NM0-6	White	Black	36"	1 Dozen
CPF-05-NS	5	•	White	•	36"	1 Dozen
BCPF-02-25-NS	2	•	White	•	36"	Bulk - 25 pack
BCPF-0230-25-NS	2	•	White	•	30"	Bulk - 25 pack
BCPF-0239-25-NS	2	•	White	•	39"	Bulk - 25 pack
BCPF-0241-25-NS	2	•	White	•	41"	Bulk - 25 pack
BCPF-02C-25-NS	2	•	White	Blue	36"	Bulk - 25 pack
BCPF-02C30-25-NS	2	•	White	Blue	30"	Bulk - 25 pack
BCPF-02C38-25-NS	2	•	White	Blue	38"	Bulk - 25 pack
BCPF-02C39-25-NS	2	•	White	Blue	39"	Bulk - 25 pack
BCPF-02C41-25-NS	2	•	White	Blue	41"	Bulk - 25 pack
BCPF-02K-25-NS	2	•	White	Black	36"	Bulk - 25 pack
BCPF-02K30-25-NS	2	•	White	Black	30"	Bulk - 25 pack
BCPF-02K38-25-NS	2	•	White	Black	38"	Bulk - 25 pack
BCPF-02K39-25-NS	2	•	White	Black	39"	Bulk - 25 pack
BCPF-02K41-25-NS	2	•	White	Black	41"	Bulk - 25 pack
BCPF-05-25-NS	5	•	White	•	36"	Bulk - 25 pack
BCPF-30-25-NS	3/0	•	White	•	36"	Bulk - 25 pack

NM0-6

26mm
1/2 Circle
Taper Point



NCCS

48mm
1/2 Circle
Conventional Cutting



Other sizes and combinations available upon request