

Grid-Loc™ & CP-HDR-Grid™

Disposable sterile template grids
for use during Prostate Brachytherapy


~Pre-sterilized, single use disposable template for patient safety.

~Eliminates cleaning time and cost, and potential for contamination from debris remaining in holes.

~Fully reimbursable
CPT Code: 77332


~ No waiting for reprocessing.

~ Competitively priced.



Grid-Loc™
(Patient Side)

CP-HDR-Grid™
(Doctor Side)



Polyurethane Coated
to facilitate catheter grip



CP Medical
Portland, Oregon
(800) 950 - 2763

Product	Color	Product Number	Size
CP-HDR-Grid	Black	CPG-015-HDR	15 Gauge Needle
Grid-Loc	Blue	CPG-015-GL-HDR	15 Gauge Needle