

MONOMID®

Monofilament Nylon
non-absorbable suture

CHARACTERISTICS

Nylon exhibits:

- Uniform diameter of the suture
- Maximum flexibility of the strand
- Excellent and consistent knot security
- Minimal memory

INDICATIONS

Monomid® suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological procedures.

COLOR

Black or Blue

TENSILE STRENGTH

As with any nylon suture, though not absorbed, progressive hydrolysis of the suture in vivo may result in gradual loss of tensile strength over time. Nylon suture should not be used where permanent retention of tensile strength is required.

PACKAGING

Monomid® is available in USP sizes 6-0 through 2 (metric sizes 0.7-5).

Monomid® is available with needles in one dozen boxes, or in veterinary suture cassettes. Three dozen suture boxes available on customer request.



803 NE 25th Avenue, Portland, OR 97232 USA

ph. 800 950-2763 • fax 503 230-9993 • www.cpmedical.com

