

POLYPRO®

Monofilament Polypropylene
non-absorbable suture

CHARACTERISTICS

Polypropylene exhibits

- Uniform diameter of the suture
- Maximum flexibility of the strand
- Excellent security with snug and flattened knots
- Minimal memory and reduced fraying
- Consistent knot strength

INDICATIONS

Polypro® polypropylene suture is indicated for use in general soft tissue approximation and/or ligation, including use in cardiovascular, ophthalmic, and neurological surgery.

COLOR

Blue

TENSILE STRENGTH

Polypro® polypropylene sutures are nonabsorbable with no significant change in strength retention is known to occur in vivo.

PACKAGING

Polypro® is available in USP sizes 8-0 through 1 (metric sizes 0.4-4).

Polypro® is available with needles in one dozen boxes. Three dozen suture boxes available on customer request.



803 NE 25th Avenue, Portland, OR 97232 USA
ph. 800 950-2763 • fax 503 230-9993 • www.cpmedical.com

